


Measuring your Dog's Neck for a Wag Bandana



Use either a soft tape or a piece of string to measure your pet's neck. Wrap the tape or string loosely around the preferred area where the bandana will sit - the thickest part of the neck - and remember not to pull it tight, just snug enough so that you can fit two fingers between it and your dog's neck.

If you're using string, just cut it at the desired length and measure it against a ruler.

If your pup hates that idea, just remove their collar, lay it flat and measure using only one side of the buckle.

Then using the table below, select your appropriate bandana size.

Bandana Size Guide			
Size	Bandana Length (approximate)	Neck Size	Example Breeds
Extra Small	51cm	18cm	Toy Dog Breeds: Chihuahuas, Toy Poodles, Maltese, or similar.
Small/Medium	54cm	27cm	Dachshunds, Yorkies, Westies, Jack Russells, Pugs, Boston Terriers and similar sized dogs.
Medium/Large	74cm	37cm	Springers, Cockers, Labradors, Boxers, Staffies and similar sized dogs.
Extra Large	98cm	47cm	Mastiffs, Great Danes, Iris Wolfhounds, Dogue De Bordeaux, Leon Berger and similar sized.